

# Haemodynamic Response to Change in Fluid Status in and Experimental Rats - Case Study



The measure of life.

## Case Study

### Haemodynamic Response To Change In Fluid Status In An Experimental Rat

Dr G Gunawandena<sup>1</sup>, Dr K Ranasinghe<sup>1</sup>, C Aitken<sup>2</sup>, R Phillips<sup>2,3</sup>

1. Therapeutic Research Unit, Princess Alexandra Hospital, Brisbane, Australia, 2. USCOM Ltd, Sydney, Australia,  
3. Department of Medicine, University of Queensland, Brisbane, Australia

#### Introduction:

An improved understanding of physiology and pathophysiology in animals and humans can be acquired by studying animal models. Cardiac Output (CO) is a core physiologic observation and essential to most animal studies, if only to confirm that an intervention or therapy is not cardiopressive. Currently, as in humans, most animal CO studies involve an invasive thermodilution method or an inaccurate alternative. USCOM is an accurate noninvasive CW Doppler which both increases sensitivity to haemodynamic change and is ethically more desirable. USCOM has been used in humans, dogs, sheep, pigs (large and small), but not previously in rats as heart rates are often above 300bpm.

#### Presentation:

1 rat was examined as part of an IRB approved study to evaluate therapy in burns. USCOM was used to evaluate haemodynamics and determine the feasibility of USCOM measurement of high HR circulation.

#### Observations:

Baseline rat weight: 455gms

#### Protocol

1. Baseline post anaesthesia (n=10)
2. Post insertion of tracheal tube and invasive lines (n=5)
3. Pre-burn (n=10)
4. 1 minute post burn (n=10)
5. 5 minutes post burn. Hypovolemia from fluid shift, and low BP (n=11)
6. 0.6ml of fluid infused (n=10)
7. 1 minute post fluid infusion (n=43)
8. SV Variability study (n=5)
9. Post fluid CO measures (n=10)



**Uscom Limited**  
ABN 35 091 028 090  
Level 7, 10 Loftus Street  
Sydney NSW 2000 Australia  
T +612 9247 4144 F +612 9247 8157  
[www.uscom.com.au](http://www.uscom.com.au)

# Haemodynamic Response to Change in Fluid Status in and Experimental Rats - Case Study



The measure of life.




Figure 1 USCOM screen during rat study demonstrating an average HR of 402bpm

## USCOM values:

- Baseline post anaesthesia (n=10)  
SV=0.62±0.05mls    HR=329±26    CO=0.21±0.02l/min
- Post insertion of tracheal tube and invasive lines (n=5)  
SV=0.62±0.05mls    HR=358±23    CO=0.22±0.02l/min
- Pre-burn (n=10)  
SV=0.62±0.12mls    HR=339±31    CO=0.21±0.03l/min
- 1 minute post burn (n=10)  
SV=0.47±0.03mls    HR=402±30    CO=0.19±0.02l/min
- 5 minutes post burn. Hypovolemia from fluid shift, low BP, (n=11)  
SV=0.45±0.07mls    HR=357±27    CO=0.16±0.03l/min
- 0.6ml of fluid infused (n=10)  
SV=0.59±0.11mls    HR=343±38    CO=0.20±0.02l/min
- 1 minute post fluid infusion (n=43)  
SV=0.63±0.13mls    HR=375±36    CO=0.23±0.05l/min
- SV Variability study (n=5)  
SV=0.72±0.07mls    HR=358±18    CO=0.26±0.02l/min
- Post fluid CO measures (n=10)  
SV=0.69±0.13mls    HR=356±47    CO=0.24±0.04l/min
- Total Observations (n=114)  
SV=0.60±0.13mls



# Haemodynamic Response to Change in Fluid Status in and Experimental Rats - Case Study



The measure of life.




Figure 2 SV, HR and CO with CIs at each intervention over full time course



# Haemodynamic Response to Change in Fluid Status in and Experimental Rats - Case Study



**Uscom**

*The measure of life.*




Figure 3 Scatterplot of all SV measures over full time course

	r	p	n
HR vs SV	-0.346	0.000	114
HR vs CO	0.165	0.080	114
SV vs CO	0.863	0.000	114

Table 1 Correlation of HR with SV and CO, demonstrating HR is a poor correlate of CO and SV

## Observations:

HR is not a good surrogate of SV or CO.

SV rose with fluid infusion from 0.63cm<sup>3</sup> to 0.72cm<sup>3</sup> (14%) and then fell during equilibration to 0.69cm<sup>3</sup> (10%).

SV changes associated with fluid depletion and infusion were detected using USCOM. Haemodynamic measurements were feasible at HRs up to 461bpm.



**Uscom Limited**  
 ABN 35 091 028 090  
 Level 7, 10 Loftus Street  
 Sydney NSW 2000 Australia  
 T +612 9247 4144 F +612 9247 8157  
 www.uscom.com.au

# Haemodynamic Response to Change in Fluid Status in and Experimental Rats - Case Study



*The measure of life.*

## Conclusion:

Measurement of cardiac performance in a rat was feasible using USCOM at heart rates up to 461bpm.

USCOM measured anticipated SV, HR and CO changes associated with induced hypovolemia and fluid therapy.

Beat to beat haemodynamics provides unique physiologic insights.

## References:

1. Phillips RA, Hood SG, Jacobson BM, Lichtenthal PR, Burstow DJ, West MJ, May CN. Measurement of CO by flow probe, USCOM and PAC in conscious sheep at rest and after dobutamine. (Abstract) ANZICS ASM, October 21st 2005
2. Lester A Critchley, Zhi Y Peng, Benny S Fok, Jules Flach, Simon C Wong, Anna Lee, Robert A Phillips. Testing the reliability of a new Ultrasonic Cardiac Output Monitor, the USCOM, using aortic flow probes in anaesthetized dogs. *Anesth Analg* 2005;100:748-53



**Uscom Limited**  
ABN 35 091 028 090  
Level 7, 10 Loftus Street  
Sydney NSW 2000 Australia  
T +612 9247 4144 F +612 9247 8157  
[www.uscom.com.au](http://www.uscom.com.au)