

## VALIDATION OF USCOM CO MEASUREMENTS IN PRETERM NEONATES BY COMPARISON WITH ECHOCARDIOGRAPHY

Phillips RA, Paradisis M, Evans NJ, Southwell DL, Burstow DJ, West MJ. The University of Queensland, Brisbane, Australia, Royal Prince Alfred Hospital, Sydney, Australia, and The Prince Charles Hospital, Brisbane, Australia.

---

---

**INTRODUCTION:** Measurement of cardiac output (CO) in pre-term neonates is important for clinical management. Doppler ultrasound is the preferred method for measurement of CO however the pulmonary and aortic diameters for calculating flow volumes are small, and measurement using 2D ultrasound requires expertise and experience, particularly of the pulmonary valve. The USCOM (USCOM Ltd, Sydney, Australia) is a novel 2D independent device which uses CW Doppler and anthropometrics to determine flow volumes. The device is simpler to operate and less expensive than conventional echo.

This study was to compare 2D echo and USCOM CO measurements in pre-term neonates.

**METHOD:** After IRB approval 66 paired measures of transpulmonary CO were acquired in 37 pre-term neonates (weight  $1.13 \pm 0.47$ kg) using conventional echocardiography, combining 2D and CW Doppler, and the USCOM device. Signals were acquired and analysed independently and in a blinded fashion, and values compared by two tailed t-tests and Bland-Altman bias analysis.

**RESULTS:** Mean values of transpulmonary CO were  $0.36 \pm 0.19$ l/min by echo and  $0.37 \pm 0.14$ l/min by USCOM and not significantly different ( $p < 0.005$ ). The mean difference between measures was  $0.00 \pm 0.08$ l/min, with a mean % error of -3.7%. The smaller SD associated with USCOM convert to smaller 95% CIs and a possible increased sensitivity for detection of haemodynamic change.

**CONCLUSION:** These results suggest that USCOM is as accurate for measurement of neonatal CO as conventional echo, and may make a cost-effective contribution for neonatal haemodynamic management.